

Staebel Slideezi Support Straps - Pair

staebel 

Product Code: 210064




An aid to reduce the risk of over-stretching and reaching in frequent supine transfers.

These Support Straps are used to reduce the risk of over-stretching and reaching in supine transfers. No need to launder reducing the chance of losses.

A new concept in support straps, these can be wiped down and/or washed in soapy water, rinsed and then left to dry. Used to reduce the risk of over-stretching and reaching in frequent supine transfers, they simply clip onto the handles of the flat Slideezi with handles and other patient handling equipment.

- Sold in packs of two so practitioners can work in pairs to position patients safely and reduce the risk of carer injury
- Durable strap construction means the straps can take the weight of most patients
- Easy-to-use buckle mechanism allows carers to attach the straps to various lateral transfer boards
- Easy to wipe down



Supplied
in a pair



Easy to
Clean



Durable
Design



Simple
to Use

Specifications

Product Code	210064
Width	50 mm
Length	750 mm
SWL	100 kg

